

体外诊断原料
首选合作伙伴

Dedication to the IVD industry?

We I.V.DO that™

Medix Biochemica 2023 技术讲座讲义

讲座时间：2023 年 5 月 29 日 10:00 - 12:00

讲座地点：南昌绿地国际博览中心一楼 A2 会议室



应您所需，量身定制 ——分子诊断用工程化DNA聚合酶 Shaping DNA polymerases for your needs

RAMON KRANASTER
博士/研发总监

2023年5月29日 南昌
Medix Biochemica




目录 Contents

1. 天然酶及其优化方式
Nature's enzymes and how they can be improved
2. 酶的定向进化策略
Description of Directed evolution strategy
3. 工程化酶的应用实例
Examples for applications with improved enzymes
4. 如何让您的物流更安全高效
One more thing: How to relax your logistics team

DNA聚合酶是生物行业基础技术的关键角色

DNA polymerases are key-players in fundamental techniques in Biotechnology




测序、检测、基因分型、克隆等。
Sequencing, Detection, Genotyping, Cloning, etc.

Medix Biochemica

DNA聚合酶 — 天然的“分子机器”

DNA polymerases — nature's molecular machines

- 延伸速度 Velocity:
35-150 nts/sec
- 错配率 Error rates:
 1.1×10^{-4} errors/bp - 8.9×10^{-5} errors/bp
- 碱基切除修复 (BER) 和核苷酸切除修复 (NER) 的校准和错配修复, 总错配率大约是 10^{-8} - 10^{-10} errors/bp。
Proofreading and mismatch repair by base and nucleotide excision repair (BER and NER) overall error-rates of 10^{-8} - 10^{-10} errors/bp.
- 在真核细胞中, 错配率约 $>10^{-10}$ 。
In eukaryotic cells error-rates $>10^{-10}$ are reached.
- 人类每个细胞基因组中含有大约60亿个核苷酸!
In living, human cells: 6 billion nucleotides!



Medix Biochemica

DNA聚合酶 — 天然的“分子机器”

DNA polymerases — nature's molecular machines

“大自然并未完全根据我们的需求塑造它”
“Nature didn't shape them to our needs”

我们的贡献：
Our contribution:


- 为生物技术应用量身定制DNA聚合酶
- Tailoring DNA polymerases for biotec applications
- 具有新特性或改良特性的DNA聚合酶突变体
- DNA polymerase mutant with new or improved features
- 使您的生物技术应用更简单、更快、更可靠
- easier, faster, more reliable biotechnological application

5

Medix Biochemica

我们的贡献：为生物技术应用量身定制DNA聚合酶

Our contribution: Tailoring DNA polymerases for biotec applications



6


Medix Biochemica

方法 Approach

通过诱导随机突变...
by random mutagenesis...

Medix Biochemica

随机突变 ——DNA聚合酶的工程化改造 Engineering of DNA polymerases by random mutagenesis



Medix Biochemica

Kranaster, R. and Marx, A. (2010), Engineered DNA Polymerases in Biotechnology. ChemBioChem, 11: 2077–2084. doi: 10.1002/cbic.201000215

Gloeckner, C., Kranaster, R. and Marx, A. (2010), Directed Evolution of DNA Polymerases: Construction and Screening of DNA Polymerase Mutant Libraries. Current Protocols in Chemical Biology, 2: 89–109.


方法 Approach

通过人工设计…
by rational designs…

Medix Biochemica

9

人工设计
——DNA聚合酶的工程化改造
Engineering of DNA polymerases by
rational designs



10

Medix Biochemica

HiDi® DNA polymerase

改造实例 I Example I

高鉴别能力
高灵敏度
等位基因特异性PCR的金标准

用于突变检测的特化功能酶
A specialized enzyme for mutation detections


Medix Biochemica

11

改造实例 I — HiDi® DNA polymerase Example I — HiDi® DNA polymerase

HiDi® 意味着“High Discrimination”，即高鉴别能力。
HiDi® is standing for “High Discrimination”.

Primer design




12


Medix Biochemica

改造实例 I — 通过 HiDi® 实现高度选择性的等位基因特异性PCR Example I — Highly selective allele-specific PCR with HiDi®


Primer design
rs72921001
3' -ACAAGTATCTTACCCCTATCTCTCCCTAGT TTT T GACATCC -5'
C-allele primer: 5' -GAATGGGATAGAGAAGGGATCAAAAAG-3'
A-allele primer: 5' -GAATGGGATAGAGAAGGGATCAAAAAT-3'
Reverse primer: 5' -CTGCTGCTGAAAATGGATTGTG-3'




Taq DNA polymerase
(N公司)



Platinum Taq DNA polymerase
(TF公司)



HiDi® DNA polymerase
(myPOLS Biotec)



HiDi® DNA polymerase 是本实验中唯一能够准确区分引物/模板错配的酶。
HiDi® DNA polymerase is the only enzyme which separates right from wrong primer/template.

Medix Biochemica


13

改造实例 I — HiDi® 是液体活检应用的理想选择 Example I — HiDi® is ideal for liquid biopsy approaches


Mutation detection test: BRAF c.1799T>A (V600E), rs113488022, homo sapiens

target sequence: 5' . . . ATAGGTGATTTTGGTCTAGCTACAGT / A GAAATCTCGATG . . .
forward primer: 5' -GGTGATTTTGGTCTAGCTACAGA-3'
reverse primer: 5' -ACCATCCACAAAATGGATCCA-3'
taqman-probe: 5' -ROX-TCGATGGAGTGGGTCCCATCAGTTT-BMNQ590-3'


Taq DNA Polymerase
Hot Start Version (TKR公司)



GoTaq® hot start
DNA Polymerase (P公司)



HiDi Taq DNA polymerase
(myPOLS Biotec, #9201)



HiDi® Taq DNA polymerase 能够在 >10,000 个野生型序列中检测到单个突变拷贝。
HiDi® Taq DNA polymerase enables detection of single mutations copies in the presence of >10.000 wildtype sequences.

Medix Biochemica

14

Volcano® RT-PCR Master Mix

改造实例 II Example II

一种具有PCR反应活性、
耐高温的逆转录酶

A PCR active, thermostable reverse transcriptase

Medix Biochemica


15

改造实例 II — 兼具DNA聚合酶活性、耐高温的逆转录酶？

Example II — a PCR active, thermostable reverse transcriptase?


传统的
一步法 RT-qPCR

conventional
one-step RT-qPCR protocol



零步法 RT-qPCR

0-step RT-qPCR protocol



我们可以使RT-PCR更简单、更快吗？
Can we make RT-PCR easier and shorter?


16

Medix Biochemica



17

Volcano[®] DNA polymerase ...



RNA -> DNA


... 95°C 半衰期 >40 min

... has a half-life at 95°C of >40 min

18

真正热稳定性 — 兼具DNA聚合酶活性、耐高温的逆转录酶

True thermostability — a PCR-active thermostable reverse transcriptase



在更高的温度下进行逆转录有利于处理复杂样品 (如基质、RNA二级结构), 并提高特异性。
 Performing reverse transcription at a higher temperature can benefit assays dealing with complex samples (e.g. matrices, RNA secondary structures) and increase specificity



19

Blatter, N., Bergen, K., Nolte, O., Welte, W., Diederichs, K., Mayer, J., Wieland, M., and Marx, A. (2013). Structure and Function of an RNA-Reading Thermostable DNA Polymerase. *Angew. Chem. Int. Ed.*, 52, 11935–11939. doi: 10.1002/anie.201306655
 R. Kranaster, M. Drum, N. Engel, M. Weidmann, F. T. Hufert and A. Marx. (2010). One-step RNA pathogen detection with reverse transcriptase activity of a mutated thermostable Thermus aquaticus DNA polymerase. *Biotechnol. J.*, 5(2), 224-31.
 Kranaster R, Zeller J, Kühn B, Marx A. (2016) Neues Enzym mit Reverse Transkriptase- und DNA-Polymerase-Funktion, *BioSpektrum*, 2, 164-165.

Medix Biochemica

零步法RT-qPCR — 使用 Volcano® RT-PCR Master Mix 实例

0-step RT-qPCR — examples with Volcano® RT-PCR Master Mix



Volcano2G RT-PCR mastermix with GreenDye
 qRT-PCR-Assay, hGAPDH mRNA
 total RNA dilution series: 100-0.01 ng/reaction
 Sybr Green I channel

hGAPDH-Assay
 hGAPDH primer fwd 5'-GAAGGTGAAGGTCGGAGTCAAC-3'
 hGAPDH primer rev 5'-GCTTCCCCTTCTCAGCCTTG-3'

RT-PCR-Protocol: 95°C 60 sec
 95°C 15 sec (40 cycles)
 60°C 60 sec

melting: 95°C 10 sec
 65°C 60 sec
 97°C 1 sec ramp 0.2°C/sec, 5 readings/°C

- 无需单独的反转录步骤
- simplified RT-PCR without isothermal RT-steps
- 也适用于染料法
- dsDNA dye compatible
- 定量结果更可靠
- reliable quantification

20


Medix Biochemica

零步法RT-qPCR — 使用 Volcano® RT-PCR Master Mix 实例

0-step RT-qPCR — examples with Volcano® RT-PCR Master Mix

Volcano2G RT-PCR mastermix
 IDT PrimeTime qPCR Assay Hs.PT.58.25447000 (from Integrated DNA Technologies),
 NONO mRNA, FAM-Channel
 total RNA dilution series: 100-0.01 ng/reaction

IDT PrimeTime qPCR Assay Hs.PT.58.25447000, NONO mRNA
 NONO primer fwd 5'-TGATGCTTTCTTGGAGTATGGT-3'
 NONO primer rev 5'-AGAGGAGAAGTCGAGGTTAGAG-3'
 NONO probe 5'-/6-
 FAM/ACTGGGAGG/ZEN/CACTTTGCTGTCTG/IABkFQ/-3'
 RT-PCR-Protocol: 95°C 1 min
 95°C 10 sec
 60°C 1 min (50 cycles)



➔ 适用于探针法 (Taqman)
 compatible with Taqman probes

21


Medix Biochemica

零步法RT-qPCR — 细胞样本直扩 (无需抽提!)

0-step RT-qPCR — direct from Cell samples (no clean-up!)

细胞mRNA直扩:
mRNA detection direct from cells:
 Taqman probe-based RT-qPCR, HPRT1 mRNA, Cy5-Channel
 fresh HEK-293 dilution series in water: 5000 – 10 cells/reaction,
 (cell were detached, resuspended and diluted according to the manual)

HPRT1- Assay
 HPRT1 primer fwd 5'-AGACTTTGCTTTCCTTGGTCAG-3'
 HPRT1 primer rev 5'-TCAAGGGCATATCCTACAACAA-3'
 HPRT1 probe 5'-Cy5-AAGCTTGCTGGTGAAGGA-BHQ2-3'
 RT-PCR-Protocol: 95°C 1 min
 95°C 15 sec
 60°C 1 min (40 cycles)



➔ 无需抽提——细胞样本直扩 PCR
 It works also without extraction – directPCR from Cells

22

Medix Biochemica

Isotherm3G DNA polymerase

兼具逆转录酶和DNA聚合酶活性
高特异性
快速检测

改造实例 III Example III

兼具逆转录活性的 LAMP DNA 聚合酶
A LAMP-capable DNA polymerase with reverse transcription activity


Medix Biochemica

23

兼具反转录活性的工程化 LAMP 酶

An engineered LAMP enzyme with improved RT-activity

等温扩增
Isothermal amplification





➔ 同 Bst DNA 聚合酶一样，可检测 DNA 模板并得到可靠的检测结果
DNA targets are reliably detected as with any other Bst DNA polymerase

24

Medix Biochemica

兼具反转录活性的工程化 LAMP 酶


An engineered LAMP enzyme with improved RT-activity



无需添加逆转录酶，可检测 RNA 模板、并在**更短时间内**获得可靠的检测结果。
RNA targets are reliably detected in **shorter time** with our engineered polymerase
No extra addition of a reverse transcriptase is needed!

兼具反转录活性的工程化 LAMP 酶

An engineered LAMP enzyme with improved RT-activity



无需添加逆转录酶，可检测 RNA 模板并获得可靠的检测结果。
与其他供应商的酶 (聚合酶3.0版) 相比，无非特异性扩增。
RNA targets are reliably detected with our engineered polymerase
A competitor product ("Polymerase 3.0") yields unspecific, unreliable signals

PlexTaq®
5x qPCR Multiplex Master Mix

粗样本直扩
高耐抑制性
多重检测

改造实例 IV Example IV

高耐抑制性、可直接用于血液样本qPCR的
预混液


Inhibitor resistant polymerases a PCR mastermix
for real-time DirectPCRs from blood samples

Medix Biochemica

27

传统的基因分型方法较单一

Standard Genotyping approach is tedious



28

Medix Biochemica

传统的基因分型方法较单一

Standard Genotyping approach is tedious



- 节省成本, 无需提取试剂或设备
No costs for extraction reagents and instruments
- 更低的出错率
Less error-prone
- 更高的通量
Higher throughput

Medix Biochemica

29

血液直扩、实时定量PCR — PlexTaq® 5x qPCR Multiplex Master Mix Real-time DirectPCRs from Blood samples – PlexTaq® 5x qPCR Multiplex

- **FAM 通道:** 人类凝血风险因子5-SNP基因靶点
FAM channel: Human coagulation risk factor 5- SNP gene target
- **HEX 通道:** spike蛋白对照 (10^4 c/rxn)
HEX channel: spike IC control target (10^4 c/rxn)
- **Cy5 通道:** RNaseP靶基因 (人血液样本中存在的) - 靶基因随着血液样本量的变化而自然变化
Cy5 channel: RNaseP target (human gene present in human blood sample) – target will naturally vary with amount of blood sample
- 不同血液样本的稀释 (EDTA、柠檬酸盐和肝素处理) 0.1%-20% 终浓度
Dilution of different blood samples (EDTA, citrate and Heparin treated)
From 0.1% - 20% final concentration

PCR protocol:

Denaturation and lysis of samples	95 °C	2 min (is meant for denaturing and lysis of samples)	
Denaturation	95 °C	10 sec	} 50x
Annealing	60 °C	40 sec	
(Fluorescence measurement in FAM, HEX and Cy5 channel)			
Melting curve measurement in FAM channel: 45 - 85 °C			






我们的预混液可以省略耗时的提取步骤, 进行PCR直扩。
Our PCR mix allows directPCRs without time-consuming extraction steps

Medix Biochemica

30




PlexTaq® 5x qPCR Multiplex Master Mix, 血液样本直扩, 未受任何抑制影响
PlexTaq® 5x qPCR Multiplex doesn't show any inhibition in blood samples



- HEX通道: spike蛋白对照 (10⁴ c/rxn)
 HEX channel: spike IC control target (10⁴ c/rxn)
- Taqman 探针
 Taqman probe

Medix Biochemica

使用 PlexTaq® 进行基因分型PCR, 血液样本直扩, 无需抽提
Blood samples extraction steps are not needed for Genotyping PCRs with PlexTaq®




- 杂交探针
 Hybridisation probes
- FAM通道:
 FAM channel:
 人类凝血风险因子5-SNP基因靶点
 Human coagulation risk factor 5- SNP gene target

Medix Biochemica

PlexTaq[®] 预混液用于多重病原体检测

Multiplex pathogen detections with- PlexTaq[®] 5x qPCR Multiplex



➔ PlexTaq[®] 可以支持高度多重检测应用，结果稳定均一、无明显偏差。
PlexTaq[®] allows high multiplexing and uniform amplification without a strong PCR bias

Medix Biochemica

33

Taq Hotstart DNA Polymerase

室温下高度稳定
高耐抑制性
多重检测

改造实例 V
Example V

可于室温下进行反应配置
——别担心，我们已搞定！

Engineered DNA polymerase for room-temperature reaction setups – don't worry we have it covered!

Medix Biochemica

34

将包含引物/探针的酶预混液，于室温下孵育

Incubation of a full mastermix including primers and probes at room-temperature

其他供应商Taq DNA 聚合酶经21°C过夜孵育：
信号强度减弱、灵敏度降低



Competitors, standard
Taq DNA polymerase:

Loss of signal intensities and reduced limit of
detections

工程化热启动 Taq DNA 聚合酶经21°C过夜孵育：
信号强度及灵敏度未受任何影响

Our engineered,
Taq Hotstart DNA Polymerase:

Don't worry anymore about room-temperature setups!



更多的可能...

One more thing...

运输不畅时，你会担心产品的质量吗？

Your shipment is stuck, delayed or returned?

每个冻干微球可包含多对引物/探针
Every bead can contain multiple primers and probes

室温运输或储存
Shipping and storage at room-temperature

仅需添加样品即可
Only sample needs to be added



LyoBeads 冻干微球
(1 Bead = 1 rxn) (1冻干微球=1个反应)



让我们的冻干技术为您赋能！

Make use of our freeze-drying capabilities!

冻干粉原料用于多次反应：
1份冻干粉=25次反应
Multiple reactions LyoCakes:
For example 1 Cake = 25 rxns

联系我们获得定制化服务！

Establishment and production of customized products possible. Get in touch!





总结 Summary

为生物技术应用量身定制DNA聚合酶
Tailoring DNA polymerases for biotec applications.

- 随机突变——DNA聚合酶的工程化改造
Engineering of DNA polymerases by random mutagenesis
- 人工设计——DNA聚合酶的工程化改造
Engineering of DNA polymerases by rational designs


专业，应您所需——赋能分子诊断
Your needs are our expertise - make molecular diagnosis capable.

- 突变检测 (HiDi®)
Mutation detections (HiDi®)
- 耐高温RT-PCR (Volcano®)
Hot RT-PCR (Volcano®)
- 适用于DNA和RNA模板的通用LAMP酶 (Isotherm3G)
One LAMP for all (Isotherm3G)
- 血液样本直扩、无需抽提 (PlexTaq® & 血液直扩PCR)
Skip Blood extraction steps (PlexTaq® & DirectBlood)

让我们的冻干技术为您赋能
Make use of our freeze-drying capabilities

客户定制服务
Other customized service

Medix Biochemica



Thank you

Meet us at booth A6-2609

Medix Biochemica

联系我们：

地址：上海市闵行区浦江镇绿洲环路 10 号 6 幢 11 层

电话：021-68119180,68119105,68119181

邮箱：medixchina@medixbiochemica.com

Medix Biochemica

www.cnmedixbiochemica.com

