# Medix Biochemica

## Anti-h HE4 4503 SPTN-5

Product overview			
Catalog number	100608		
Specificity	Antibody recognizes human epididymis protein 4 (HE4)		
Description	Monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components.		
Product buffer solution	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN $_3$ as a preservative		
Shelf life and storage	24 months from manufacturing at 2–8 °C		
Subclass	lgG₁		
Analyte description	Human epididymis protein 4 (HE4), also known as WAP four-disulfide core domain protein 2, is a 124 amino acid long protease inhibitor. Serum HE4 is often measured together with CA125 to monitor progression of epithelial ovarian cancer after treatment.		
Parameters tested on e	each lot		
Product appearance	Liquid, may turn slightly opaque during storage		
Product concentration	5.0 mg/ml (+/- 10 %)		
Immunoreactivity	80–120 % compared to the reference sample in an FIA test		
IEF Profile	6.1–7.1		
Purity	≥ 95 %		
Kinetic parameters			
Accordiation rate constant	Nat Determined (N/D)		

Association rate constant	Not Determined (N/D)

Dissociation rate constant N/D

Affinity constant N/D

Determination antigen

**Determination method** 



medix@medixbiochemica.com www.medixbiochemica.com VAT reg. no. FI14631532

#### Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.

# Medix Biochemica

### **Cross-reactivities**

**Epitope** 

### Does not recognize cancer antigens CA125, CA19-9 or CA15-3

		DETECTION				
		4501	4502	4503	4505	4506
CAPTURE	4501	-	-	+	+	+
	4502	-	-	+	+	+
	4503	-	-	-	-	-
	4505	-	-	-	-	-
	4506	-	-	-	-	-

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Platforms tested	FIA	
Antigens tested	N/D	
Product stability	temperatures affect the antigen bindi	ife given on the first page is based on
Miscellaneous	-	

References



**Oy Medix Biochemica Ab** Klovinpellontie 3, Fl-02180 Espoo, Finland medix@medixbiochemica.com www.medixbiochemica.com VAT reg. no. FI14631532

#### Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.