



# **Product specifications**

| Name                    | Anti-hCG beta 5004 SP-1  |
|-------------------------|--|
| Specificity             | Antibody recognizes human chorionic gonadotropin and its free beta subunit   |
| Description             | Monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components.  |
| Product code            | 100004   |
| Product buffer solution | 0.9 % NaCl, 0.095 % NaN $_3$ as a preservative   |
| Shelf life and storage  | 36 months from manufacturing at 2–8 °C   |
| Subclass                | IgG <sub>1</sub>   |
| Analyte description     | Human chorionic gonadotropin (hCG) is a glycoprotein hormone produced in pregnancy by<br>the developing embryo soon after conception and later by the syncytiotrophoblast (part of<br>the placenta). Its role is to prevent the disintegration of the corpus luteum of the ovary and<br>thereby maintain progesterone production that is critical for a pregnancy in humans. Early<br>pregnancy testing, in general, is based on the detection of hCG. hCG is produced also by some<br>tumors, but it is not known whether this production is a contributing cause or an effect of<br>tumorigenesis. |

## Parameters tested on each lot

| Product appearance    | Liquid, may turn slightly opaque during storage          |
|-----------------------|--|
| Product concentration | 1.0 mg/ml (+/- 10 %)                                     |
| Immunoreactivity      | 80–120 % compared to the reference sample in an FIA test |
| IEF Profile           | 6.5–7.6  |
| Purity                | ≥ 95 %   |

### **Kinetic parameters**

| Association rate constant  | hCG: 2.0 x 10 <sup>6</sup> 1/Ms and hCG $\beta$ : 5.0 x 10 <sup>5</sup> 1/Ms   |
|----------------------------|--|
| Dissociation rate constant | hCG: 3.7 x 10 <sup>-5</sup> 1/s and hCG $\beta$ : 1.2 x 10 <sup>-4</sup> 1/s   |
| Affinity constant          | hCG: K <sub>A</sub> = 5.3 x 10 <sup>10</sup> 1/M; K <sub>D</sub> = 1.9 x 10 <sup>-11</sup> M (= 0.02 nM)<br>hCGβ: K <sub>A</sub> = 4.3 x 10 <sup>9</sup> 1/M; K <sub>D</sub> = 2.3 x 10 <sup>-10</sup> M (= 0.23 nM) |
| Determination method       | SPR analysis (ProteOn XPR36)   |
| Determination antigen      | hCG, Scripps (Cat C0714, Lot 2430801); hCGβ, Scripps (Cat C0914, Lot 2310001)  |

**Oy Medix Biochemica Ab** Klovinpellontie 3 FI-02180 Espoo, Finland

medix@medixbiochemica.com www.medixbiochemica.com VAT reg.no. FI14631532

Legal disclaimer Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. Medix/MAB and Medix/Antigens are trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.





#### **Cross-reactivities**

LH (K<sub>D</sub> = 17 nM) Does not recognize hCGa, FSH, or TSH

Epitope

Beta-2 as described in Berger et al. (2013). The antibody recognizes both intact hCG and free  $\beta$  subunit.

#### Pair recommendations

|         |                | DETECTION |      |      |      |      |                |      |               |      |      |      |
|---------|----------------|-----------|------|------|------|------|----------------|------|---------------|------|------|------|
|         |                | hCG beta  |      |      |      |      |                |      | alpha subunit |      |      |      |
|         |                | 5004      | 5006 | 5008 | 5009 | 5011 | 5012<br>free β | 5014 | 5016          | 5501 | 5503 | 6601 |
|         | 5004           | -         | -    | -    | +    | +    | -              | +    | -             | +    | +    | +    |
|         | 5006           | -         | -    | -    | -    | -    | -              | +    | -             | +    | +    | +    |
| JRE     | 5008           | -         | -    | -    | +    | -    | -              | +    | -             | +    | +    | +    |
|         | 5009           | +         | +    | +    | -    | -    | -              | +    | +             | -    | -    | +    |
| CAPTURE | 5011           | +         | +    | +    | -    | -    | -              | +    | +             | -    | -    | +    |
| 0       | 5012<br>free β | +         | +    | +    | -    | -    | -              | +    | +             | -    | -    | -    |
|         | 5014           | +         | +    | +    | +    | +    | -              | -    | +             | +    | +    | +    |
|         | 5016           | -         | -    | -    | -    | -    | -              | +    |               | +    | +    | +    |

Following pairs are especially recommended for free hCG beta assays: CLIA: 5012 (capture) – 5004 (detection) and 5012 – 5008

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

| Platforms tested  | FIA, CLIA  |  |  |  |  |  |
|-------------------|--|--|--|--|--|--|
| Antigens tested   | 5 .  | Native hCG antigens, Lee Biosolutions 189-10 and 189-11<br>Native β-hCG antigen, Lee Biosolutions 325-11 |  |  |  |  |
| Product stability | TEMPERATURE, TIME<br>Not Determinend (N/D)   | RESULT<br>-  |  |  |  |  |
| Miscellaneous     | -  |  |  |  |  |  |
| References        | Berger, P., Paus, E., Hemken, P.M., Sturgeon, C., Stewart, W.W., Skinner, J.P., Harwick, L |  |  |  |  |  |

Berger, P., Paus, E., Hemken, P.M., Sturgeon, C., Stewart, W.W., Skinner, J.P., Harwick, L.C., Saldana, S.C., Ramsay, C.S., Rupprecht, K.R., Olsen, K.H., Bidart, J.M. and Stenman, U.H. (2013) Candidate epitopes for measurement of hCG and related molecules: the second ISOBM TD-7 workshop. Tumor Biol., 34: 4033-4057.

**Oy Medix Biochemica Ab** Klovinpellontie 3 FI-02180 Espoo, Finland medix@medixbiochemica.com www.medixbiochemica.com VAT reg.no. FI14631532 Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB and MedixAntigens are trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.